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
Manufacturer	Tanita Upper Arm Type Blood Pressure Monitor
Measurement Method	Oscillometric
Measurement Range	Pressure: 0 – 300mmHg, Pulse: 40 – 199 pulses/minute
Scale	Pressure: 1mmHg
Accuracy	Pressure: ± 3 mmHg, Pulse: Read value $\pm 4\%$
Display Information	Diastolic pressure, systolic pressure and pulse, simultaneous display
Memory Function	Saves 60 measurement results (diastolic pressure, systolic pressure and pulse)
Battery Life	Approx. 350 measurements
Power	DC6V, AA alkaline dry-cell batteries $\times 4$
Pressurization	Automatic pressurization via pump
Decompression	Constant decompression through automatic exhaust valve system
Exhaust	Automatic electromagnetic valve exhaust system
Safety Mechanism	Rapid exhaust at 300mmHg via electromagnetic valve
Electric Shock Protection	Internal power supply
Worn Parts	BF type equipment
Main Unit Dimensions	D150 \times W128 \times H70 (mm)
Weight	Approx. 255g (including batteries)
Operational Temperature	Temperature +10 – +40 $^{\circ}\text{C}$ / Humidity 85% or less (w/o condensation)
Transportation/Storage Temperature	Temperature -5 – +50 $^{\circ}\text{C}$ / Humidity 85% or less (w/o condensation)
Supported Arm Size	22 – 33 cm.

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EMC Tables

Guidance and manufacturer's declaration-electromagnetic emissions

The Tanita Upper Arm Type Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Tanita Upper Arm Type Blood Pressure Monitor should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The Tanita Upper Arm Type Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any
RF emissions CISPR 11	Class B	The Tanita Upper Arm Type Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration-electromagnetic immunity

The Tanita Upper Arm Type Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Tanita Upper Arm Type Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration-electromagnetic immunity

The Tanita Upper Arm Type Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Tanita Upper Arm Type Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Tanita Upper Arm Type Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d=2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. </p>

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Tanita Upper Arm Type Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Tanita Upper Arm Type Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Tanita Upper Arm Type Blood Pressure Monitor.

EMC Tables (continued)

Recommended separation distances between portable and mobile RF communications equipment and the Tanita Upper Arm Type Blood Pressure Monitor

The Tanita Upper Arm Type Blood Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Tanita Upper Arm Type Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Tanita Upper Arm Type Blood Pressure Monitor as recommended below,

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected